



# Coronavirus Disinfectant Selector Guide

| Product Name & Number                  | Active | Dilution Ratio | pH       | Acinetobacter | Enterobacteriaceae (carbapenem-resistance) | Hepatitis A | Hepatitis B | Hepatitis C | HIV    | Human Coronavirus | Human Coronavirus (SARS Associated) | Influenza | Klebsiella | MRSA   | Norovirus | Pseudomonas aeruginosa | Staphylococcus aureus | VRE    | Tuberculosis (TB) | VISA   | VRSA   | Canine Parvo | Salmonella enterica/choleeraesuis |
|--|--------|----------------|----------|---------------|--|-------------|-------------|-------------|--------|-------------------|-------------------------------------|-----------|------------|--------|-----------|------------------------|-----------------------|--------|-------------------|--------|--------|--------------|-----------------------------------|
| Triforce (333)                         | Quat   | 1:256          | Alkaline | 3 min         |  |             | 3 min       | 3 min       | 1 min  | 3 min             | 10 min                              | 1 min     | 3 min      | 3 min  |           | 3 min                  | 3 min                 | 3 min  |                   | 3 min  |        |              | 3 min                             |
| Quat-Stat 5™ (341)                     | Quat   | 1:256          | Alkaline | 5 min         | 5 min                                      |             | 5 min       | 5 min       | 5 min  | 5 min             | 5 min                               | 5 min     | 5 min      | 5 min  | 5 min     | 5 min                  | 5 min                 | 5 min  |                   |        | 5 min  |              | 5 min                             |
| Fight-Bac® RTU (311)                   | Quat   | RTU            | Alkaline |               |  | 10 min      | 1 min       | 1 min       | 1 min  | 2 min             | 2 min                               | 2 min     |            | 3 min  | 30 sec    | 3 min                  | 3 min                 | 3 min  | 5 min             | 3 min  |        | 10 min       | 3 min                             |
| pH7Q (316)                             | Quat   | 1:64           | Neutral  | 10 min        |  |             | 10 min      | 10 min      | 4 min  | 1 min             | 10 min                              | 2 min     | 10 min     | 10 min |           | 10 min                 | 10 min                | 10 min |                   | 10 min |        |              | 10 min                            |
| pH7Q Dual (355)                        | Quat   | 1:256          | Neutral  | 10 min        | 10 min                                     |             | 10 min      | 10 min      | 10 min | 10 min            | 10 min                              | 10 min    | 10 min     | 10 min |           | 10 min                 | 10 min                | 10 min |                   | 10 min | 10 min | 10 min       | 10 min                            |
| Pine Quat (304)                        | Quat   | 1:32           | Neutral  | 10 min        |  |             | 10 min      | 10 min      | 4 min  | 10 min            | 10 min                              | 10 min    | 10 min     | 10 min |           |                        | 10 min                | 10 min |                   | 10 min |        |              | 10 min                            |
| Symplicity™ Sanibet™ Multi-Range (237) | Quat   | 3 oz/5 gal     | Neutral  |               |  |             | 10 min      | 10 min      | 10 min | 10 min            | 10 min                              | 10 min    | 10 min     | 10 min | 10 min    | 10 min                 | 10 min                |        |                   |        |        |              | 10 min                            |

- Triforce Sanitizer for soft, non-food contact surfaces at 1:256